
Roche CustomBiotech – IVT reaction starting condition

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Non-medical device/reagent, PMR_21112022_000120

IVT Reaction protocol

IVT reaction protocol (starting conditions)

The following conditions may be a good start:

- 10x Transcription Buffer (basic buffer): 400 mM Tris-HCl (pH 8.0), 20 mM spermidine, 100 mM DTT
- Reaction set up: 20 uL reaction volume; 500 ng template; 1 hr incubation at 37 C
- Reagents used (all Roche CB reagents):
 - 1 µL T7 RNA Polymerase
 - 10 mM each NTP
 - 1 µL PPase
 - 20 U RNase Inhibitor
 - 40 mM MgCl₂

Nucleotide and Mg²⁺ concentration has a high influence on the reaction. As a general guideline it is better to keep the ratio of NTP:Mg²⁺ between 0.7-1, during the titration. Alternative 10x Transcription Buffer (not Tris based) : 1 M HEPES-KOH (pH 7.5), 20 mM Spermidine, 400 mM DTT.

When optimizing the conditions, the customers may vary the ratio of the NTPs from 1:1:1:1 to something slightly different. Often, GTP is used in somewhat higher concentration.

UTP can be replaced entirely by N1-Me-Pseudo-UTP (Me-psUTP), or a certain percentage of UTP along with Me-psUTP. This really depends on careful optimization. The amount needed of each NTP depends on the target sequence and further optimizing the process to maximize mRNA yield.