

Summary

Expiry date, Material number and lot numbering system

Affected Products All Roche CustomBiotech products

Description

Expiry Date:

Roche CustomBiotech products are delivered with an expiry date in the format yyyy-mm. The shelf life is given until the last day of the specified month.

Material number:

All materials supplied by Roche Diagnostics GmbH are identified by a material number with an 11-digits format.

The first 8-digits represent the material and the last 3 digits represent the designated use.

Please note: Some labels/documents may contend only the first 8 digits of the material number.

Lot number:

Each batch of a material supplied by Roche Diagnostics GmbH can be identified by a lot number with an 8-digits format.

The first 6-digits represent the mother lot which is unique and sequentially generated by our SAP system. The last two digits provide specific information about the filling of the material.

If there are any questions or concerns with respect to this notification, please feel free to contact your local key account manager.

Best regards,

i.V.

**Dr. Heike
Koenig**

Digital signiert von Dr. Heike
Koenig
DN: cn=Dr. Heike Koenig, c=DE,
email=heike.koenig@roche.com
Datum: 2022.06.28 09:47:47
+02'00'

Heike König
Sr. Director, Head Global Product
Management CustomBiotech

i.V.

ehsess

Digital signiert von ehsses
DN: cn=ehsses,
email=sylvia.ehses@roche.com
Datum: 2022.06.28 12:50:05
+02'00'

Sylvia Ehse
Head of Portfolio Management & Supply B2B,
Diagnostics Operations Penzberg